

Antibodies Product Datasheet

Product Name:	PGC-1 α Monoclonal Antibody
Catalog Number:	N0518
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human PGC-1 α expressed in E. Coli.
Concentration:	
Gene Name:	PPARGC1A
Protein Name:	Peroxisome proliferator-activated receptor gamma coactivator 1- α p
Other Name:	PPARGC1A; LEM6; PGC1; PGC1A; PPARGC1; Peroxisome proliferator-activated receptor gamma coactivator 1- α ; PGC-1- α ; PPAR-gamma coactivator 1- α ; PPARGC-1- α ; Ligand effect modulator 6
Human Swiss Prot No.:	Q9UBK2
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	PGC-1 α Monoclonal Antibody detects endogenous levels of PGC-1 α protein.
Storage and Stability:	-20°C/1 year