

Antibodies Product Datasheet

Product Name:	Pinch-1 Monoclonal Antibody
Catalog Number:	N0519
Reactivity:	Human
Applications:	WB,IF/ICC,FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Pinch-1 expressed in E. Coli.
Concentration:	
Gene Name:	LIMS1
Protein Name:	LIM and senescent cell antigen-like-containing domain protein 1
Other Name:	LIMS1; PINCH; PINCH1; LIM and senescent cell antigen-like-containing domain protein 1; Particularly interesting new Cys-His protein 1; PINCH-1; Renal carcinoma antigen NY-REN-48
Human Swiss Prot No.:	P48059
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Pinch-1 Monoclonal Antibody detects endogenous levels of Pinch-1 protein.
Storage and Stability:	-20°C/1 year