

Antibodies Product Datasheet

Product Name:	Pirh2 Monoclonal Antibody
Catalog Number:	N0520
Reactivity:	Human,Rat
Applications:	WB,IHC-p,IF/ICCFCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Pirh2 expressed in E. Coli.
Concentration:	
Gene Name:	RCHY1
Protein Name:	RING finger and CHY zinc finger domain-containing protein 1
Other Name:	RCHY1; ARNIP; CHIMP; PIRH2; RNF199; ZNF363; RING finger and CHY zinc finger domain-containing protein 1; Androgen receptor N-terminal-interacting protein; CH-rich-interacting match with PLAG1; E3 ubiquitin-protein ligase Pirh2; RING finger
Human Swiss Prot No.:	Q96PM5
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Pirh2 Monoclonal Antibody detects endogenous levels of Pirh2 protein.
Storage and Stability:	-20°C/1 year