

Antibodies Product Datasheet

Product Name:	PLAP Monoclonal Antibody
Catalog Number:	N0521
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human PLAP expressed in E. Coli.
Concentration:	
Gene Name:	ALPP
Protein Name:	Alkaline phosphatase, placental type
Other Name:	ALPP; PLAP; Alkaline phosphatase; placental type; Alkaline phosphatase Regan isozyme; Placental alkaline phosphatase 1; PLAP-1
Human Swiss Prot No.:	P05187
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	PLAP Monoclonal Antibody detects endogenous levels of PLAP protein.
Storage and Stability:	-20°C/1 year