

Antibodies Product Datasheet

Product Name:	Podocalyxin-like 1 Monoclonal Antibody
Catalog Number:	N0525
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Podocalyxin-like 1 expressed in E. Coli.
Concentration:	
Gene Name:	PODXL
Protein Name:	Podocalyxin
Other Name:	PODXL; PCLP; PCLP1; Podocalyxin; GCTM-2 antigen; Gp200; Podocalyxin-like protein 1; PC; PCLP-1
Human Swiss Prot No.:	O00592
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Podocalyxin-like 1 Monoclonal Antibody detects endogenous levels of Podocalyxin-like 1 protein.
Storage and Stability:	-20°C/1 year