

Antibodies Product Datasheet

Product Name:	PR3 Monoclonal Antibody
Catalog Number:	N0531
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human PR3 expressed in E. Coli.
Concentration:	
Gene Name:	PRTN3
Protein Name:	Myeloblastin
Other Name:	PRTN3; MBN; Myeloblastin; AGP7; C-ANCA antigen; Leukocyte proteinase 3; PR-3; PR3; Neutrophil proteinase 4; NP-4; P29; Wegener autoantigen
Human Swiss Prot No.:	P24158
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	PR3 Monoclonal Antibody detects endogenous levels of PR3 protein.
Storage and Stability:	-20°C/1 year