

Antibodies Product Datasheet

Product Name:	PRAK Monoclonal Antibody
Catalog Number:	N0532
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of PRAK expressed in E. Coli.
Concentration:	
Gene Name:	MAPKAPK5
Protein Name:	MAP kinase-activated protein kinase 5
Other Name:	MAPKAPK5; PRAK; MAP kinase-activated protein kinase 5; MAPK-activated protein kinase 5; MAPKAP kinase 5; MAPKAP-K5; MAPKAPK-5; MK-5; MK5; p38-regulated/activated protein kinase; PRAK
Human Swiss Prot No.:	Q8IW41
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	PRAK Monoclonal Antibody detects endogenous levels of PRAK protein.
Storage and Stability:	-20°C/1 year