

Antibodies Product Datasheet

Product Name:	PRK2 Monoclonal Antibody
Catalog Number:	N0533
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IHC-p, IF (paraffin section), FCM, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human PRK2 expressed in E. Coli.
Concentration:	
Gene Name:	PKN2
Protein Name:	Serine/threonine-protein kinase N2
Other Name:	PKN2; PRK2; PRKCL2; Serine/threonine-protein kinase N2; PKN gamma; Protein kinase C-like 2; Protein-kinase C-related kinase 2
Human Swiss Prot No.:	Q16513
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	PRK2 Monoclonal Antibody detects endogenous levels of PRK2 protein.
Storage and Stability:	-20°C/1 year