

## Antibodies Product Datasheet

|                               |                                                                                                                            |
|-------------------------------|----------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | Rab 10 Monoclonal Antibody                                                                                                 |
| <b>Catalog Number:</b>        | N0543                                                                                                                      |
| <b>Reactivity:</b>            | Human, Mouse                                                                                                               |
| <b>Applications:</b>          | WB, IF/ICC, ELISA                                                                                                          |
| <b>Dilutions:</b>             | Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000.<br>ELISA: 1/10000. Not yet tested in other applications. |
| <b>Modification:</b>          | Unmodified/ Total                                                                                                          |
| <b>Source:</b>                | Monoclonal Mouse                                                                                                           |
| <b>Purification:</b>          | Affinity purification                                                                                                      |
| <b>Immunogen:</b>             | Purified recombinant fragment of human Rab 10 expressed in E. Coli.                                                        |
| <b>Concentration:</b>         |                                                                                                                            |
| <b>Gene Name:</b>             | RAB10                                                                                                                      |
| <b>Protein Name:</b>          | Ras-related protein Rab-10                                                                                                 |
| <b>Other Name:</b>            | RAB10; Ras-related protein Rab-10                                                                                          |
| <b>Human Swiss Prot No.:</b>  | P61026                                                                                                                     |
| <b>Formulation:</b>           | Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.                                                       |
| <b>Specificity:</b>           | Rab 10 Monoclonal Antibody detects endogenous levels of Rab 10 protein.                                                    |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                               |