

Antibodies Product Datasheet

Product Name:	RAG-2 Monoclonal Antibody
Catalog Number:	N0549
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human RAG-2 (350-527aa) expressed in E. Coli.
Concentration:	
Gene Name:	RAG2
Protein Name:	V(D)J recombination-activating protein 2
Other Name:	RAG2; V(D)J recombination-activating protein 2; RAG-2
Human Swiss Prot No.:	P55895
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	RAG-2 Monoclonal Antibody detects endogenous levels of RAG-2 protein.
Storage and Stability:	-20°C/1 year