

Antibodies Product Datasheet

Product Name:	Rak Monoclonal Antibody
Catalog Number:	N0550
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of Rak (aa2-300) expressed in E. Coli.
Concentration:	
Gene Name:	FRK
Protein Name:	Tyrosine-protein kinase FRK
Other Name:	FRK; PTK5; RAK; Tyrosine-protein kinase FRK; FYN-related kinase; Nuclear tyrosine protein kinase RAK; Protein-tyrosine kinase 5
Human Swiss Prot No.:	P42685
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Rak Monoclonal Antibody detects endogenous levels of Rak protein.
Storage and Stability:	-20°C/1 year