

Antibodies Product Datasheet

Product Name:	Reg Iα Monoclonal Antibody
Catalog Number:	N0554
Reactivity:	Human
Applications:	IHC-p,IF/ICC,ELISA
Dilutions:	Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Reg Iα fused with hlgGFc tag expressed in HEK293 cell line.
Concentration:	
Gene Name:	REG1A
Protein Name:	Lithostathine-1-alpha
Other Name:	REG1A; PSPS; PSPS1; REG; Lithostathine-1-alpha; Islet cells regeneration factor; ICRF; Islet of Langerhans regenerating protein; REG; Pancreatic stone protein; PSP; Pancreatic thread protein; PTP; Regenerating islet-derived protein 1-alpha;
Human Swiss Prot No.:	P05451
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Reg Iα Monoclonal Antibody detects endogenous levels of Reg Iα protein.
Storage and Stability:	-20°C/1 year