

Antibodies Product Datasheet

Product Name:	Ret Monoclonal Antibody
Catalog Number:	N0555
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of Ret (aa896-1063) expressed in E. Coli.
Concentration:	
Gene Name:	RET
Protein Name:	Proto-oncogene tyrosine-protein kinase receptor Ret
Other Name:	RET; CDHF12; CDHR16; PTC; RET51; Proto-oncogene tyrosine-protein kinase receptor Ret; Cadherin family member 12; Proto-oncogene c-Ret
Human Swiss Prot No.:	P07949
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Ret Monoclonal Antibody detects endogenous levels of Ret protein.
Storage and Stability:	-20°C/1 year