

Antibodies Product Datasheet

Product Name:	Ribosomal Protein S27 Monoclonal Antibody
Catalog Number:	N0557
Reactivity:	Human
Applications:	IHC-p,IF(paraffin section),ELISA
Dilutions:	Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of Ribosomal Protein S27 expressed in E. Coli.
Concentration:	
Gene Name:	RPS27
Protein Name:	40S ribosomal protein S27
Other Name:	RPS27; MPS1; 40S ribosomal protein S27; Metallopan-stimulin 1; MPS-1
Human Swiss Prot No.:	P42677
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Ribosomal Protein S27 Monoclonal Antibody detects endogenous levels of Ribosomal Protein S27 protein.
Storage and Stability:	-20°C/1 year