

Antibodies Product Datasheet

Product Name:	Ron Monoclonal Antibody
Catalog Number:	N0561
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Ron (aa210-320) expressed in E. Coli.
Concentration:	
Gene Name:	MST1R
Protein Name:	Macrophage-stimulating protein receptor
Other Name:	MST1R; PTK8; RON; Macrophage-stimulating protein receptor; MSP receptor; CDw136; Protein-tyrosine kinase 8; p185-Ron; CD antigen CD136
Human Swiss Prot No.:	Q04912
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Ron Monoclonal Antibody detects endogenous levels of Ron protein.
Storage and Stability:	-20°C/1 year