

Antibodies Product Datasheet

Product Name:	ROR1 Monoclonal Antibody
Catalog Number:	N0562
Reactivity:	Human
Applications:	WB,IF/ICC,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Recombinant extracellular fragment of human ROR1 (aa30-406) fused with hlgGFc tag, expressed in HEK293 cells
Concentration:	
Gene Name:	ROR1
Protein Name:	Tyrosine-protein kinase transmembrane receptor ROR1
Other Name:	ROR1; NTRKR1; Tyrosine-protein kinase transmembrane receptor ROR1; Neurotrophic tyrosine kinase; receptor-related 1
Human Swiss Prot No.:	Q01973
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	ROR1 Monoclonal Antibody detects endogenous levels of ROR1 protein.
Storage and Stability:	-20°C/1 year