

Antibodies Product Datasheet

Product Name:	Rsk-2 Monoclonal Antibody
Catalog Number:	N0564
Reactivity:	Human
Applications:	WB,IF/ICC,FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Rsk-2 expressed in E. Coli.
Concentration:	
Gene Name:	RPS6KA3
Protein Name:	Ribosomal protein S6 kinase alpha-3
Other Name:	RPS6KA3; ISPK1; MAPKAPK1B; RSK2; Ribosomal protein S6 kinase alpha-3; S6K-alpha-3; 90 kDa ribosomal protein S6 kinase 3; p90-RSK 3; p90RSK3; Insulin-stimulated protein kinase 1; ISPK-1; MAP kinase-activated protein kinase 1b; MAPK-activated
Human Swiss Prot No.:	P51812
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Rsk-2 Monoclonal Antibody detects endogenous levels of Rsk-2 protein.
Storage and Stability:	-20°C/1 year