

Antibodies Product Datasheet

Product Name:	Rtn-3 Monoclonal Antibody
Catalog Number:	N0566
Reactivity:	Human
Applications:	IHC-p,IF(paraffin section),ELISA
Dilutions:	Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of Rtn-3 expressed in E. Coli.
Concentration:	
Gene Name:	RTN3
Protein Name:	Reticulon-3
Other Name:	RTN3; ASYIP; NSPL2; Reticulon-3; Neuroendocrine-specific protein-like 2; NSP-like protein 2; Neuroendocrine-specific protein-like II; NSP-like protein II; NSPLII
Human Swiss Prot No.:	O95197
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Rtn-3 Monoclonal Antibody detects endogenous levels of Rtn-3 protein.
Storage and Stability:	-20°C/1 year