

Antibodies Product Datasheet

Product Name:	Skp2 p45 Monoclonal Antibody
Catalog Number:	N0578
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of Skp2 p45 (aa1-130) expressed in E. Coli.
Concentration:	
Gene Name:	SKP2
Protein Name:	S-phase kinase-associated protein 2
Other Name:	SKP2; FBXL1; S-phase kinase-associated protein 2; Cyclin-A/CDK2-associated protein p45; F-box protein Skp2; F-box/LRR-repeat protein 1; p45skp2
Human Swiss Prot No.:	Q13309
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Skp2 p45 Monoclonal Antibody detects endogenous levels of Skp2 p45 protein.
Storage and Stability:	-20°C/1 year