

Antibodies Product Datasheet

Product Name:	SLUG Monoclonal Antibody
Catalog Number:	N0579
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human SLUG expressed in E. Coli.
Concentration:	
Gene Name:	SNAI2
Protein Name:	Zinc finger protein SNAI2
Other Name:	SNAI2; SLUG; SLUGH; Zinc finger protein SNAI2; Neural crest transcription factor Slug; Protein snail homolog 2
Human Swiss Prot No.:	O43623
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	SLUG Monoclonal Antibody detects endogenous levels of SLUG protein.
Storage and Stability:	-20°C/1 year