

Antibodies Product Datasheet

Product Name:	Smad6 Monoclonal Antibody
Catalog Number:	N0584
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Smad6 expressed in E. Coli.
Concentration:	
Gene Name:	SMAD6
Protein Name:	Mothers against decapentaplegic homolog 6
Other Name:	SMAD6; MADH6; Mothers against decapentaplegic homolog 6; MAD homolog 6; Mothers against DPP homolog 6; SMAD family member 6; SMAD 6; Smad6; hSMAD6
Human Swiss Prot No.:	O43541
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Smad6 Monoclonal Antibody detects endogenous levels of Smad6 protein.
Storage and Stability:	-20°C/1 year