

Antibodies Product Datasheet

Product Name:	SNAI 1 Monoclonal Antibody
Catalog Number:	N0588
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human SNAI 1 expressed in E. Coli.
Concentration:	
Gene Name:	SNAI1
Protein Name:	Zinc finger protein SNAI1(snail)
Other Name:	SNAI1; SNAH; Zinc finger protein SNAI1; Protein snail homolog 1; Protein sna
Human Swiss Prot No.:	O95863
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	SNAI 1 Monoclonal Antibody detects endogenous levels of SNAI 1 protein.
Storage and Stability:	-20°C/1 year