

Antibodies Product Datasheet

Product Name:	Synuclein- α Monoclonal Antibody
Catalog Number:	N0605
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of Synuclein- α expressed in E. Coli.
Concentration:	
Gene Name:	SNCA
Protein Name:	Alpha-synuclein
Other Name:	SNCA; NACP; PARK1; Alpha-synuclein; Non-A beta component of AD amyloid; Non-A4 component of amyloid precursor; NACP
Human Swiss Prot No.:	P37840
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Synuclein- α Monoclonal Antibody detects endogenous levels of Synuclein- α protein.
Storage and Stability:	-20°C/1 year