

Antibodies Product Datasheet

Product Name:	Synuclein- γ Monoclonal Antibody
Catalog Number:	N0606
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of Synuclein- γ expressed in E. Coli.
Concentration:	1 mg/ml
Gene Name:	SNCG
Protein Name:	Gamma-synuclein
Other Name:	SNCG; BCSG1; PERSYN; PRSN; Gamma-synuclein; Breast cancer-specific gene 1 protein; Persyn; Synoretin; SR
Human Swiss Prot No.:	O76070
Formulation:	Purified antibody in PBS containing 0.03% sodium azide.
Specificity:	Synuclein- γ Monoclonal Antibody detects endogenous levels of Synuclein- γ protein.
Storage and Stability:	-20°C/1 year