

Antibodies Product Datasheet

Product Name:	TBL1X Monoclonal Antibody
Catalog Number:	N0609
Reactivity:	Human,Monkey
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human TBL1X expressed in E. Coli.
Concentration:	
Gene Name:	TBL1X
Protein Name:	F-box-like/WD repeat-containing protein TBL1X
Other Name:	TBL1X; TBL1; F-box-like/WD repeat-containing protein TBL1X; SMAP55; Transducin beta-like protein 1X; Transducin-beta-like protein 1; X-linked
Human Swiss Prot No.:	O60907
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	TBL1X Monoclonal Antibody detects endogenous levels of TBL1X protein.
Storage and Stability:	-20°C/1 year