

Antibodies Product Datasheet

Product Name:	TCP-1 β Monoclonal Antibody
Catalog Number:	N0612
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IF/ICC, FCM, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human TCP-1 β expressed in E. Coli.
Concentration:	
Gene Name:	CCT2
Protein Name:	T-complex protein 1 subunit beta
Other Name:	CCT2; 99D8.1; CCTB; T-complex protein 1 subunit beta; TCP-1-beta; CCT-beta
Human Swiss Prot No.:	P78371
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	TCP-1 β Monoclonal Antibody detects endogenous levels of TCP-1 β protein.
Storage and Stability:	-20°C/1 year