

Antibodies Product Datasheet

Product Name:	TH-POK Monoclonal Antibody
Catalog Number:	N0618
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human TH-POK expressed in E. Coli.
Concentration:	
Gene Name:	ZBTB7B
Protein Name:	Zinc finger and BTB domain-containing protein 7B
Other Name:	ZBTB7B; ZBTB15; ZFP67; ZNF857B; Zinc finger and BTB domain-containing protein 7B; Krueppel-related zinc finger protein cKrox; hckKrox; T-helper-inducing POZ/Krueppel-like factor; Zinc finger and BTB domain-containing protein 15; Zinc finger
Human Swiss Prot No.:	O15156
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	TH-POK Monoclonal Antibody detects endogenous levels of TH-POK protein.
Storage and Stability:	-20°C/1 year