

Antibodies Product Datasheet

Product Name:	TNK1 Monoclonal Antibody
Catalog Number:	N0622
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of TNK1 (aa451-560) expressed in E. Coli.
Concentration:	
Gene Name:	TNK1
Protein Name:	Non-receptor tyrosine-protein kinase TNK1
Other Name:	TNK1; Non-receptor tyrosine-protein kinase TNK1; CD38 negative kinase 1
Human Swiss Prot No.:	Q13470
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	TNK1 Monoclonal Antibody detects endogenous levels of TNK1 protein.
Storage and Stability:	-20°C/1 year