

Antibodies Product Datasheet

Product Name:	TORC2 Monoclonal Antibody
Catalog Number:	N0624
Reactivity:	Human, Monkey
Applications:	WB, IHC-p, IF/ICCFCM, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human TORC2 expressed in E. Coli.
Concentration:	
Gene Name:	CRTC2
Protein Name:	CREB-regulated transcription coactivator 2
Other Name:	CRTC2; TORC2; CREB-regulated transcription coactivator 2; Transducer of regulated cAMP response element-binding protein 2; TORC-2; Transducer of CREB protein 2
Human Swiss Prot No.:	Q53ET0
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	TORC2 Monoclonal Antibody detects endogenous levels of TORC2 protein.
Storage and Stability:	-20°C/1 year