

Antibodies Product Datasheet

Product Name:	TORC3 Monoclonal Antibody
Catalog Number:	N0625
Reactivity:	Human, Monkey
Applications:	WB, IHC-p, IF/ICELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human TORC3 expressed in E. Coli.
Concentration:	
Gene Name:	CRTC3
Protein Name:	CREB-regulated transcription coactivator 3
Other Name:	CRTC3; TORC3; CREB-regulated transcription coactivator 3; Transducer of regulated cAMP response element-binding protein 3; TORC-3; Transducer of CREB protein 3
Human Swiss Prot No.:	Q6UUV7
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	TORC3 Monoclonal Antibody detects endogenous levels of TORC3 protein.
Storage and Stability:	-20°C/1 year