

Antibodies Product Datasheet

Product Name:	Trk A Monoclonal Antibody
Catalog Number:	N0627
Reactivity:	Human
Applications:	WB,IF/ICC,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant extracellular fragment of human Trk A (aa33-423) fused with hIgGFc tag expressed in HEK293 cell line.
Concentration:	
Gene Name:	NTRK1
Protein Name:	High affinity nerve growth factor receptor
Other Name:	NTRK1; MTC; TRK; TRKA; High affinity nerve growth factor receptor; Neurotrophic tyrosine kinase receptor type 1; TRK1-transforming tyrosine kinase protein; Tropomyosin-related kinase A; Tyrosine kinase receptor; Tyrosine kinase receptor A;
Human Swiss Prot No.:	P04629
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Trk A Monoclonal Antibody detects endogenous levels of Trk A protein.
Storage and Stability:	-20°C/1 year