

Antibodies Product Datasheet

Product Name:	TudorSN Monoclonal Antibody
Catalog Number:	N0634
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of TudorSN (aa361-485) expressed in E. Coli.
Concentration:	
Gene Name:	SND1
Protein Name:	Staphylococcal nuclease domain-containing protein 1
Other Name:	SND1; TDRD11; Staphylococcal nuclease domain-containing protein 1; 100 kDa coactivator; EBNA2 coactivator p100; Tudor domain-containing protein 11; p100 co-activator
Human Swiss Prot No.:	Q7KZF4
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	TudorSN Monoclonal Antibody detects endogenous levels of TudorSN protein.
Storage and Stability:	-20°C/1 year