

Antibodies Product Datasheet

Product Name:	TUG Monoclonal Antibody
Catalog Number:	N0635
Reactivity:	Mouse
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of TUG expressed in E. Coli.
Concentration:	1 mg/ml
Gene Name:	ASPSCR1
Protein Name:	Tether containing UBX domain for GLUT4
Other Name:	ASPSCR1; ASPL; RCC17; TUG; UBXD9; UBXN9; Tether containing UBX domain for GLUT4; Alveolar soft part sarcoma chromosomal region candidate gene 1 protein; Alveolar soft part sarcoma locus; Renal papillary cell carcinoma protein 17; UBX domain
Human Swiss Prot No.:	
Formulation:	Purified antibody in PBS containing 0.03% sodium azide.
Specificity:	TUG Monoclonal Antibody detects endogenous levels of TUG protein.
Storage and Stability:	-20°C/1 year