

Antibodies Product Datasheet

Product Name:	UBC9 Monoclonal Antibody
Catalog Number:	N0641
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCFCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human UBC9 expressed in E. Coli.
Concentration:	
Gene Name:	UBE2I
Protein Name:	SUMO-conjugating enzyme UBC9
Other Name:	UBE2I; UBC9; UBCE9; SUMO-conjugating enzyme UBC9; SUMO-protein ligase; Ubiquitin carrier protein 9; Ubiquitin carrier protein I; Ubiquitin-conjugating enzyme E2 I; Ubiquitin-protein ligase I; p18
Human Swiss Prot No.:	P63279
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	UBC9 Monoclonal Antibody detects endogenous levels of UBC9 protein.
Storage and Stability:	-20°C/1 year