

Antibodies Product Datasheet

Product Name:	WDR5 Monoclonal Antibody
Catalog Number:	N0646
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human WDR5 expressed in E. Coli.
Concentration:	
Gene Name:	WDR5
Protein Name:	WD repeat-containing protein 5
Other Name:	WDR5; BIG3; WD repeat-containing protein 5; BMP2-induced 3-kb gene protein
Human Swiss Prot No.:	P61964
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	WDR5 Monoclonal Antibody detects endogenous levels of WDR5 protein.
Storage and Stability:	-20°C/1 year