
Antibodies Product Datasheet

Product Name:	ACE1 Monoclonal Antibody
Catalog Number:	N0655
Reactivity:	Human, Mouse, Rat, Bovine
Applications:	WB
Dilutions:	Western Blot: 1/1000 - 1/2000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant human ACE1 (N-terminus) protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	ACE
Protein Name:	Angiotensin-converting enzyme
Other Name:	ACE; DCP; DCP1; Angiotensin-converting enzyme; ACE; Dipeptidyl carboxypeptidase I; Kininase II; CD antigen CD143
Human Swiss Prot No.:	P12821
Formulation:	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Specificity:	ACE1 Monoclonal Antibody detects endogenous levels of ACE1 protein.
Storage and Stability:	-20°C/1 year