
Antibodies Product Datasheet

Product Name:	ALB Monoclonal Antibody
Catalog Number:	N0659
Reactivity:	Human
Applications:	WB
Dilutions:	Western Blot: 1/1000 - 1/2000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant human ALB protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	ALB
Protein Name:	Serum albumin
Other Name:	ALB; GIG20; GIG42; Serum albumin
Human Swiss Prot No.:	P02768
Formulation:	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Specificity:	ALB Monoclonal Antibody detects endogenous levels of ALB protein.
Storage and Stability:	-20°C/1 year