

Antibodies Product Datasheet

Product Name:	ALDH3A2 Monoclonal Antibody
Catalog Number:	N0660
Reactivity:	Human, Mouse, Rat, Rabbit
Applications:	WB
Dilutions:	Western Blot: 1/1000 - 1/2000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant human ALDH3A2 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	ALDH3A2
Protein Name:	Fatty aldehyde dehydrogenase
Other Name:	ALDH3A2; ALDH10; FALDH; Fatty aldehyde dehydrogenase; Aldehyde dehydrogenase 10; Aldehyde dehydrogenase family 3 member A2; Microsomal aldehyde dehydrogenase
Human Swiss Prot No.:	P51648
Formulation:	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Specificity:	ALDH3A2 Monoclonal Antibody detects endogenous levels of ALDH3A2 protein.
Storage and Stability:	-20°C/1 year