

Antibodies Product Datasheet

Product Name:	B23 Monoclonal Antibody
Catalog Number:	N0667
Reactivity:	Human, Mouse, Rat, Bovine, Pig, Rabbit
Applications:	WB, IF/ICC
Dilutions:	Western Blot: 1/1000 - 1/2000. Immunofluorescence: 1/100 - 1/500. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant human B23 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	NPM1
Protein Name:	Nucleophosmin
Other Name:	NPM1; NPM; Nucleophosmin; NPM; Nucleolar phosphoprotein B23; Nucleolar protein NO38; Numatrin
Human Swiss Prot No.:	P06748
Formulation:	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Specificity:	B23 Monoclonal Antibody detects endogenous levels of B23 protein.
Storage and Stability:	-20°C/1 year