

Antibodies Product Datasheet

Product Name:	CHRD1 Monoclonal Antibody
Catalog Number:	N0675
Reactivity:	Human,Rat
Applications:	WB
Dilutions:	Western Blot: 1/1000 - 1/2000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant human CHRD1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	CHRD1
Protein Name:	Chordin-like protein 1
Other Name:	CHRD1; NRLN1; Chordin-like protein 1; Neuralin-1; Neurogenesin-1; Ventroptin
Human Swiss Prot No.:	Q9BU40
Formulation:	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Specificity:	CHRD1 Monoclonal Antibody detects endogenous levels of CHRD1 protein.
Storage and Stability:	-20°C/1 year