

Antibodies Product Datasheet

Product Name:	Dynactin 1 Monoclonal Antibody
Catalog Number:	N0684
Reactivity:	Human, Mouse, Bovine, Pig
Applications:	WB, IF/ICC
Dilutions:	Western Blot: 1/1000 - 1/2000. Immunofluorescence: 1/100 - 1/500. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant human Dynactin 1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	DCTN1
Protein Name:	Dynactin subunit 1
Other Name:	DCTN1; Dynactin subunit 1; 150 kDa dynein-associated polypeptide; DAP-150; DP-150; p135; p150-glued
Human Swiss Prot No.:	Q14203
Formulation:	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Specificity:	Dynactin 1 Monoclonal Antibody detects endogenous levels of Dynactin 1 protein.
Storage and Stability:	-20°C/1 year