

Antibodies Product Datasheet

Product Name:	Fyb Monoclonal Antibody
Catalog Number:	N0690
Reactivity:	Human,Rat
Applications:	WB
Dilutions:	Western Blot: 1/1000 - 1/2000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant human Fyb protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	FYB
Protein Name:	FYN-binding protein
Other Name:	FYB; SLAP130; FYN-binding protein; Adhesion and degranulation promoting adaptor protein; ADAP; FYB-120/130; p120/p130; FYN-T-binding protein; SLAP-130; SLP-76-associated phosphoprotein
Human Swiss Prot No.:	O15117
Formulation:	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Specificity:	Fyb Monoclonal Antibody detects endogenous levels of Fyb protein.
Storage and Stability:	-20°C/1 year