

Antibodies Product Datasheet

Product Name:	HAUSP Monoclonal Antibody
Catalog Number:	N0692
Reactivity:	Human, Mouse, Rat, Dog, Pig
Applications:	WB, IF/ICC
Dilutions:	Western Blot: 1/1000 - 1/2000. Immunofluorescence: 1/100 - 1/500. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant human HAUSP (C-terminus) protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	USP7
Protein Name:	Ubiquitin carboxyl-terminal hydrolase 7
Other Name:	USP7; HAUSP; Ubiquitin carboxyl-terminal hydrolase 7; Deubiquitinating enzyme 7; Herpesvirus-associated ubiquitin-specific protease; Ubiquitin thioesterase 7; Ubiquitin-specific-processing protease 7
Human Swiss Prot No.:	Q93009
Formulation:	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Specificity:	HAUSP Monoclonal Antibody detects endogenous levels of HAUSP protein.
Storage and Stability:	-20°C/1 year