

Antibodies Product Datasheet

Product Name:	HDAC3 Monoclonal Antibody
Catalog Number:	N0693
Reactivity:	Human, Mouse, Rat, Dog, Pig, Rabbit
Applications:	WB
Dilutions:	Western Blot: 1/1000 - 1/2000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant human HDAC3 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	HDAC3
Protein Name:	Histone deacetylase 3
Other Name:	HDAC3; Histone deacetylase 3; HD3; RPD3-2; SMAP45
Human Swiss Prot No.:	O15379
Formulation:	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Specificity:	HDAC3 Monoclonal Antibody detects endogenous levels of HDAC3 protein.
Storage and Stability:	-20°C/1 year