

Antibodies Product Datasheet

Product Name:	Hip Monoclonal Antibody
Catalog Number:	N0695
Reactivity:	Human,Bovine,Dog,Pig
Applications:	WB,IHC-p,IF/ICC
Dilutions:	Western Blot: 1/1000 - 1/2000. Immunohistochemistry: 1/500 - 1/1000. Immunofluorescence: 1/100 - 1/500. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant human Hip protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	ST13
Protein Name:	Hsc70-interacting protein
Other Name:	ST13; AAG2; FAM10A1; HIP; SNC6; Hsc70-interacting protein; Hip; Aging-associated protein 2; Progesterone receptor-associated p48 protein; Protein FAM10A1; Putative tumor suppressor ST13; Renal carcinoma antigen NY-REN-33; Suppression of tum
Human Swiss Prot No.:	P50502
Formulation:	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Specificity:	Hip Monoclonal Antibody detects endogenous levels of Hip protein.
Storage and Stability:	-20°C/1 year