

Antibodies Product Datasheet

Product Name:	hnRNP U Monoclonal Antibody
Catalog Number:	N0697
Reactivity:	Human,Bovine,Pig,Rabbit
Applications:	WB
Dilutions:	Western Blot: 1/1000 - 1/2000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant human hnRNP U protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	HNRNPU
Protein Name:	Heterogeneous nuclear ribonucleoprotein U
Other Name:	HNRNPU; HNRPU; SAFA; U21.1; Heterogeneous nuclear ribonucleoprotein U; hnRNP U; Scaffold attachment factor A; SAF-A; p120; pp120
Human Swiss Prot No.:	Q00839
Formulation:	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Specificity:	hnRNP U Monoclonal Antibody detects endogenous levels of hnRNP U protein.
Storage and Stability:	-20°C/1 year