

Antibodies Product Datasheet

Product Name:	HXK II Monoclonal Antibody
Catalog Number:	N0703
Reactivity:	Human, Mouse, Rat, Pig
Applications:	WB
Dilutions:	Western Blot: 1/1000 - 1/2000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant human HXK II (N-terminus) protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	HK2
Protein Name:	Hexokinase-2
Other Name:	HK2; Hexokinase-2; Hexokinase type II; HK II; Muscle form hexokinase
Human Swiss Prot No.:	P52789
Formulation:	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Specificity:	HXK II Monoclonal Antibody detects endogenous levels of HXK II protein.
Storage and Stability:	-20°C/1 year