

Antibodies Product Datasheet

Product Name:	Ikaros Monoclonal Antibody
Catalog Number:	N0705
Reactivity:	Human, Mouse, Dog
Applications:	WB
Dilutions:	Western Blot: 1/1000 - 1/2000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant human Ikaros (C-terminus) protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	IKZF1
Protein Name:	DNA-binding protein Ikaros
Other Name:	IKZF1; IK1; IKAROS; LYF1; ZNFN1A1; DNA-binding protein Ikaros; Ikaros family zinc finger protein 1; Lymphoid transcription factor LyF-1
Human Swiss Prot No.:	Q13422
Formulation:	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Specificity:	Ikaros Monoclonal Antibody detects endogenous levels of Ikaros protein.
Storage and Stability:	-20°C/1 year