

Antibodies Product Datasheet

Product Name:	LRRFIP1 Monoclonal Antibody
Catalog Number:	N0710
Reactivity:	Human
Applications:	WB
Dilutions:	Western Blot: 1/1000 - 1/2000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant human LRRFIP1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	LRRFIP1
Protein Name:	Leucine-rich repeat flightless-interacting protein 1
Other Name:	LRRFIP1; GCF2; TRIP; Leucine-rich repeat flightless-interacting protein 1; LRR FLII-interacting protein 1; GC-binding factor 2; TAR RNA-interacting protein
Human Swiss Prot No.:	Q32MZ4
Formulation:	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Specificity:	LRRFIP1 Monoclonal Antibody detects endogenous levels of LRRFIP1 protein.
Storage and Stability:	-20°C/1 year