

Antibodies Product Datasheet

Product Name:	Mi2- α Monoclonal Antibody
Catalog Number:	N0715
Reactivity:	Human, Mouse, Rat, Bovine, Dog, Pig
Applications:	WB, IHC-p, IF/ICCFCM
Dilutions:	Western Blot: 1/1000 - 1/2000. Immunohistochemistry: 1/500 - 1/1000. Immunofluorescence: 1/100 - 1/500. Flow cytometry: 1/100 - 1/200. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant human Mi2- α (C-terminus) protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	CHD3
Protein Name:	Chromodomain-helicase-DNA-binding protein 3
Other Name:	CHD3; Chromodomain-helicase-DNA-binding protein 3; CHD-3; ATP-dependent helicase CHD3; Mi-2 autoantigen 240 kDa protein; Mi2-alpha; Zinc finger helicase; hZFH
Human Swiss Prot No.:	Q12873
Formulation:	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Specificity:	Mi2- α Monoclonal Antibody detects endogenous levels of Mi2- α protein.
Storage and Stability:	-20°C/1 year